

# inviCELL™ Platelet lysate

## Packing specification

Product number: A6802P-500, A6802P-100, A6802P-50

Specifications: 500ml, 100ml, 50ml

## Storage conditions

Stored at -20°C, valid for 36 months

## Description

inviCELL™ culture supplement is a non-xenogeneic, animal serum-free, media supplement for replacing FBS (fetal bovine serum) to support cell expansion from research through clinical trials to commercial use. inviCELL™ contains abundant growth factors and cytokines necessary for research or industrial cell growth and proliferation of multiple cell types (e.g. MSCs).

## Intended use

For human ex-vivo tissue and cell culture processing applications

## Important information

Insoluble particles may form in thawed inviCELL cell culture supplement. Published research has shown that particles will not alter the performance of the product.

## Safety information

- Follow the handling instructions outlined in the Material Safety Data Sheets (MSDSs). Wear appropriate protective eyewear, clothing, and gloves.
- inviCELL™, is a cell culture supplement derived from human single donor platelets collected from healthy donors at FDA-licensed centers. Each donor has been tested using FDA-licensed tests and found nonreactive for HBsAg, Hepatitis B core antibody(anti-HBc), HIV antibody (anti-HIV-1/2), Hepatitis C antibody (anti-HCV), HTLV-1/2 antibody (anti-HTLV-1/2), Trypanosoma cruzi antibody (anti-T. cruzi), HIV-1, HCV, HBV, WNV nucleic acid testing and Syphilis microhemagglutination test. Handle in accordance with established bio-safety practices.

## MSC culture conditions

### Media:

Complete medium is comprised of a basal media (e.g.α-MEM or other supportive media) and inviCELL™

**Culture type:** Adhesion

**Culture vessels:** Cell culture plates, T-flasks, G-Rex flasks or cell culture bags

**Temperature range:** 36°C to 38°C

**Incubator atmosphere:** Humidified atmosphere of 4–6% CO<sub>2</sub>. Ensure that proper gas exchange is achieved in culture vessels.

## Precipitation in Cell Culture

- Insoluble particles may form in thawed inviCELL™, it is recommended to remove particles by centrifuge at 3,400 xg for 3~5 minutes.



- Filtering the completed medium (e.g. 5%), after inviCELL™ is diluted in the basal medium, will not affect inviCELL™ supplemented cell culture performance. However, 0.22 µm filtering is NOT recommended for 100% concentrate inviCELL™, as this may reduce 5% inviCELL™ cell culture performance.
- Repeated freeze-thaw cycles should be avoided as they may cause an increase in insoluble particles and resulting potential decrease in inviCELL™ performance.

### Protocol

- inviCELL™ shows optimal growth of MSC at 5% (v/v) in typical cell culture media, i.e. α-MEM, which contains 2mM L-Glutamine as final concentrate.
- We recommend seeding MSCs at approximately  $3 \times 10^3 \sim 6 \times 10^3$  per cm<sup>2</sup>. • For inviCELL™ has been fibrinogen-depleted and does not require the addition of heparin in the cell culture media.

### Storage

inviCELL™ product is most stable when stored frozen until needed. The recommended storage temperature is -20°C or -80°C. Thaw frozen inviCELL™ product in 37°C water bath before use. Once inviCELL™ product is thawed, it is recommended to fully use for completed medium preparation (e.g. 5%) the same day, or to divide it into single-use aliquots and store unused aliquots at -20°C or -80°C.

### Cell Lines

Bone marrow mesenchymal stem cells

Adipose tissue derived mesenchymal stem cells

Umbilical cord derived mesenchymal stem cells

Other mesenchymal stem cells

**It can only be used for scientific research. It is forbidden to use it for human, animal or other purposes.**

